

# To determine the concentration of $\text{Co}^{2+}$ ion in water sample by spectrophotometer



**Spectrophotometry.** Many substances dissolve to give coloured solutions. The higher the concentration of solute, the more light is absorbed and the less light is transmitted through the sample. Spectrophotometry is a simple technique used to measure absorbance of solutions.


**About:** When the solution contains a *specific molecule* that absorbs light at a *specific wavelength*.

The intensity of light hitting the detector after passing through a “blank” solution is measured – this is a solution that is identical to the sample but doesn’t contain the solute. **This measurement is called “ $I_0$ ”.**

Then the intensity of light hitting the detector after passing through the sample is measured. **This measurement is “ $I$ ”.**

## What is a blank in spectrophotometry?

A blank is a sample containing everything except for the significance analyte. The blank is a sample of the solvent itself.






Figure 3: Absorbance of two different compounds

Isobestic point:




Figure 4: An example of isobestic point

## What are transmittance and absorbance?

Consider monochromatic light transmitted through a solution; with an incident intensity of  $I_0$  and a transmitted intensity of  $I$ .

The transmittance,  $T$ , of the solution is defined as the ratio of the transmitted intensity,  $I$ , over the incident intensity,  $I_0$ :

$$T = \frac{I}{I_0}$$

and takes values between 0 and 1. However, it is more commonly expressed as a percentage transmittance:

$$T(\%) = 100 \frac{I}{I_0}$$

Experimental measurements are usually made in terms of **transmittance** (T), which is defined as:

$$T = I / I_0$$


where I is the **light intensity** after it passes through the sample and  $I_0$  is the **initial light intensity**. The relation between A and T is: or  
The absorbance, A, of the solution is related to the transmittance through the following relations:

$$A = \log_{10} \frac{I_0}{I}$$
$$A = -\log_{10} T$$

or

$$A = -\log T = -\log (I / I_0).$$

*Absorption of light by a sample*



The absorbance has a logarithmic relationship to the transmittance; with an absorbance of 0 corresponding to a transmittance of 100% and an absorbance of 1 corresponding to 10% transmittance.

## The Absorbance of a Solution

For each wavelength of light passing through the spectrometer, the intensity of the light passing through the reference cell is measured. This is usually referred to as  $I_0$  - that's I for Intensity.

The intensity of the light passing through the sample cell is also measured for that wavelength - given the symbol,  $I$ . If  $I$  is less than  $I_0$ , then the sample has absorbed

some of the light (neglecting reflection of light off the cuvette surface). A simple bit of math is then done in the computer to convert this into something called the absorbance of the sample - given the symbol,  $A$ . The absorbance of a transition depends on two external **assumptions**.

1. The absorbance is directly proportional to the concentration ( $c$ ) of the solution of the sample used in the experiment. (**Beer's Law**)
2. The absorbance is directly proportional to the length of the light path ( $l$ ), which is equal to the width of the cuvette. (**Lambert's law**)

**Assumption one** relates the absorbance to concentration and can be expressed as

$$A \propto c$$

The absorbance ( $A$ ) is defined via the incident intensity  $I_0$  and transmitted intensity  $I$  by

$$A = \log_{10} \left( \frac{I_0}{I} \right)$$

**Assumption two** can be expressed as

$$A \propto l$$

Combining Equations 1 and 3:

$$A \propto cl$$

This proportionality can be converted into an equality by including a proportionality constant ( $\epsilon$ ).

$$A = \epsilon cl$$

This formula is the common form of the *Beer-Lambert Law*, although it can be also written in terms of intensities:

$$A = \log_{10} \left( \frac{I_0}{I} \right) = \epsilon lc$$

The constant  $\epsilon$  is called **molar absorptivity** or **molar extinction coefficient** and is a measure of the probability of the electronic transition. On most of the diagrams you will come across, the absorbance ranges from 0 to 1, but it can go higher than that. An absorbance of 0 at some wavelength means that no light of that particular wavelength has been absorbed. The intensities of the sample and reference beam are both the same, so the ratio  $I_0/I$  is 1 and the  $\log_{10}$  of 1 is zero.

## The Beer-Lambert Law

You will find that various different symbols are given for some of the terms in the equation - particularly for the concentration and the solution length.

$$\log_{10} \frac{I_0}{I} = \epsilon l c$$

Diagram illustrating the Beer-Lambert Law equation:  $\log_{10} \frac{I_0}{I} = \epsilon l c$ . The symbols are defined as follows:

- $\epsilon$ : Greek letter, epsilon
- $l$ : length of solution the light passes through (cm)
- $c$ : concentration of solution ( $\text{mol dm}^{-3}$ )

The Greek letter epsilon in these equations is called the **molar absorptivity** - or sometimes the molar absorption coefficient. The larger the molar absorptivity, the more probable the electronic transition. In uv spectroscopy, the concentration of the sample solution is measured in  $\text{mol L}^{-1}$  and the length of the light path in cm. Thus, given that absorbance is unitless, the units of molar absorptivity are  $\text{L mol}^{-1} \text{cm}^{-1}$ .

## Limitations of the Beer-Lambert law

The linearity of the Beer-Lambert law is limited by chemical and instrumental factors. Causes of nonlinearity include:

- deviations in absorptivity coefficients at **high concentrations** ( $>0.01M$ ) due to electrostatic interactions between molecules in close proximity
- **scattering of light due to particulates** in the sample
- **fluorescence or phosphorescence** of the sample
- changes in refractive index at high analyte concentration
- shifts in chemical equilibria as a function of concentration
- non-monochromatic radiation, deviations can be minimized by using a relatively flat part of the absorption spectrum such as the maximum of an absorption band
- stray light




Figure 3 (a): Absorption spectra of Rhodamine B solutions with different concentrations in water (b) Calibration curve of Rhodamine B in the water at measured at  $\lambda_{\max}$ .

**AIM:** To determine the concentration of  $\text{Co}^{2+}$  ion in water sample by spectrophotometer.

**INSTRUMENT:** Spectrophotometer, Filter, Cobalt Nitrate

#### **PROCEDURE:**

- (1) Preparation of 0.2M standard cobalt nitrate solution: Prepare a stock solution of 0.2M cobalt nitrate (5.82 g in 100 mL distilled water). Dilute this stock solution quantitatively as to get desired concentration of 0.04M, 0.08M, 0.12M, 0.16M solution in 100 mL SMF.
- (2) Strictly follow the instructions given on the instrument.
- (3) Measure the one of the above  $\text{Co}^{2+}$  solution against water as a 'Blank' or 'Reference' using each of the filters in turn given with the instrument. The filter giving the maximum absorbance (510 nm) will be most suitable for measuring the absorbance (A) (Known as  $\lambda_{\max}$ ).
- (4) Determine the absorbance (A) (from %Transmittance) of a test solution against water as a blank. Carry out replicate measurements so as to get the corresponding readings.





- (5) Repeat the experiment for unknown  $\text{Co}^{2+}$  concentration.
- (6) Plot a graph absorbance (A) against concentration of  $\text{Co}^{2+}$  ion. A linear curve must be obtained if the Beer- Lambert's law is valid. Find the concentration of the unknown  $\text{Co}^{2+}$  solution from this curve from its absorbance.
- (7) Tabulate the observations as follows :

OBSERVATION TABLE:

| No. | Concentration<br>"V" (M) | % Transmittance<br>(T%) | Absorbance (A)<br>(From Table) |
|-----|--------------------------|-------------------------|--------------------------------|
| 1   | 0.04                     | 70 %                    | 0.18                           |
| 2   | 0.08                     | 46 %                    | 0.36                           |
| 3   | 0.12                     | 36 %                    | 0.48                           |
| 4   | 0.16                     | 22 %                    | 0.75                           |
| 5   | 0.20                     | 14 %                    | 0.93                           |
| 6   | Unknown                  | 58 %                    | 0.25                           |

GRAPH:





**Result:**

1. The straight-line nature of concentration vs absorbance at  $\lambda_{\text{max}}$  510 nm which confirms the Beer-Lambert's Law.
2. From the graph the concentration of unknown  $\text{Co}^{2+}$  solution = ..... ??? M