



**Conductometric titration
to determine the
strength of strong acid
by strong base**

Conductometric titration

- Titration in which conductance measurement are made use of in determining the end-point
- Method used to determine end point in acid-base reactions, precipitation reactions or complexation reactions



Theory

- Ohm's Law:

voltage drop (volts = J/C)
C = Coulomb

$$V = IR$$

current (amps = C/s)

resistance (Ohms)

- Resistivity, ρ :

-- a material property that is independent of sample size and geometry

$$\rho = \frac{RA}{l}$$

surface area
of current flow

current flow
path length


- Conductivity, σ

$$\sigma = \frac{1}{\rho}$$


Depends on no. of free ions, charge & size of the ions, mobility of the ions

Titration curve of a strong acid (HCl) with a strong base (NaOH)

- Before NaOH is added, the conductance is high due to the presence of highly mobile H^+
- When the NaOH is added, the conductance falls due to the replacement of H^+ by the added Na^+ as H^+ ions react with OH^- ions to form un-dissociated water



- This decrease in the conductance continues till the equivalence point
- At the equivalence point, the solution contains only NaCl
- After the equivalence point, the conductance increases due to the large conductivity of OH^- ions
- You will get a titration curve - point of intersection of the two lines gives the point of neutralization i.e. equivalence point



AIM: Conductometric titration to determine the strength of strong acid (HCl) by strong base (NaOH)

Procedure:


- Rinse the conductivity cell with double distilled water
- Pipette out 20 mL HCl (unknown conc.) in a beaker and dip the conductivity cell in it
- Add small amount of NaOH (0.1N) from burette, stir and measure the conductance
- Measure conductance after each addition
- Take at least five readings after the end point



Observation Table

Vol. of NaOH added (V_2 , mL)	Conductance (S)

Plot Graph



Let, V_2 be the volume of 0.1N (N_2) NaOH required to complete the neutralization (end point volume from the graph) 20 mL (V_1) of x N (N_1) HCl

Calculations:


$$N_1 V_1 = N_2 V_2 \quad , \quad N_1 = N_2 V_2 / V_1$$

$$\begin{aligned} \text{Strength(g/L)} &= \text{Normality} \times \text{Eq.wt.} \\ &= N_1 \times 36.5 \end{aligned}$$

For example:

OBSERVATION TABLE:


No.	Volume of NaOH added (V ₂) (ml)	Conductance (S)
1	0	3.34×10^{-3}
2	1	3.15×10^{-3}
3	2	2.69×10^{-3}
4	3	2.43×10^{-3}
5	4	2.07×10^{-3}
6	5	1.87×10^{-3}
7	6	1.57×10^{-3}
8	7	1.29×10^{-3}
9	8	1.08×10^{-3}
10	9	0.84×10^{-3}
11	10	0.91×10^{-3}
12	11	1.05×10^{-3}
13	12	1.15×10^{-3}
14	13	1.28×10^{-3}
15	14	1.37×10^{-3}



PLOT GRAPHS:-

(i) Conductance \rightarrow Volume of NaOH added (mL)

Scale + x-axis unit = 1 ml.
y-axis unit = 0.25×10^{-3}



Results:

Normality of HCl = _____ N

Strength of HCl = _____ g/L