


# Titrimetric determination of l-Ascorbic acid (Vitamin C)


# Titration method - Redox Titration Using Iodine Solution

**Iodimetry: A direct titration with only 1 reaction:**





Reducing agent

Oxidizing agent




## Theory: Iodimetric titration



None of these substances have a distinctive color.  
**How will we know how much  $I_2$  to add?**

After all the vitamin C has reacted, **excess  $I_2$**  begins to react with starch indicator solution to form a **colored complex**



**Reagents:** Freshly boiled & cooled water, 0.1N iodine solution, 1 N sulfuric acid, 1% Starch solution

**Procedure Sample:** 120 mg of Ascorbic Acid

**Titrimetric system:**

Mode : Direct titration

Titrant : 0.1 N  $I_2$  solution


Endpoint detection : Visual, using indicator (1% starch solution)

Blank solution : 50 mL of freshly boiled & cooled water, 10 mL of 1N sulfuric acid, and 3 mL of 1% starch solution

**Procedure**


**Blank estimation:** Take blank solution & titrate against 0.1 N  $I_2$  solution. End point marked by appearance of blue color. Let the volume of titrant consume be B


**Actual estimation:** Dissolve the given sample (120g of Ascorbic acid) in 50 mL of freshly boiled & cooled water. Then, add 10 mL of 1N sulfuric acid, and 3 mL of 1% starch solution. Titrate against 0.1 N  $I_2$  solution until a persistent blue color is obtained. Let the volume of titrant consume be V. Then, calculate the percentage of ascorbic acid in the given sample



## OBSERVATION :

Blank Estimation :				Actual Estimation			
Burette : 0.1 N $I_2$ Solution				Burette : 0.1 N $I_2$ Solution			
Pipette : 10ml. 1 N $H_2SO_4$ solution + 50 ml. Freshly boiled and cooled water.				Pipette : Weigh sample + 50 ml. Freshly boiled and cooled water + 10ml 1 N $H_2SO_4$			
Indicator : 1% starch solution. (1/4 TT)				Indicator : 1% starch solution. (1/4 TT)			
Endpoint : Appearance of blue color				Endpoint : Appearance of blue color			
Concurrent Burette Reading = A..... mL				Concurrent Burette Reading = B..... mL			
I.B.R.	F.B.R.	Diff	Concurrent	I.B.R.	F.B.R.	Diff	Concurrent
0	0.1	0.1		0	10.1	10.1	
0	0.1	0.1	0.1	0	10	10	10.1
0	0.1	0.1		0	10.1	10.1	
0	0.1	0.1		0	10.1	10.1	





1 eq  $\text{I}_2 \equiv$  1 eq of Ascorbic acid

### CALCULATION :

$$1\text{N} = \frac{1 \text{ gm.eq.}}{1 \text{ L}}$$

1000 ml 1 N  $\text{I}_2$  Solution = 88 gm of / - ascorbic acid

1 ml 0.1 N  $\text{I}_2$  Solution = 8.8 mg of / - ascorbic acid

(B-A) ml 0.1 N  $\text{I}_2$  Solution. = 8.8 X (B-A) mg of / - ascorbic acid

D = \_\_\_\_\_ mg of / - ascorbic acid

Now percentage purity of /-ascorbic acid =  $100 \times \text{D} / \text{wt of sample}$

### RESULT:

Percentage purity of give sample of l-Ascorbic acid =  $\frac{D}{120} \times 100$