


Iodometric Determination of Cu in Brass Sample Solution

IODINE TITRATIONS


Iodimetry: A direct titration with only 1 reaction:



Iodometry is the titration of iodine (I_2) produced when an oxidizing analyte is added to excess I^- (iodide).

Then the iodine (I_2) is usually titrated with standard **thiosulfate** solution.


Iodometry: Not a direct titration because there are 2 reactions:




Iodometric titrations

Iodine titrations in which some oxidizing agent liberate I_2 from an I^- solution and then liberated I_2 is titrated with a standard solution of reducing agent


Estimation of $CuSO_4$, $K_2Cr_2O_7$, $KMnO_4$, Fe^{3+} , H_2O_2 , Br_2 , Cl_2 etc.



Titration Setup



- **Detection of the end point – using chemical Indicator**
- In this titration a solution of **starch** is used as indicator
- Starch reacts with iodine in the presence of iodide to form an intensely blue-colored complex, which is visible at very low concentrations of iodine



Aim: Estimation of Cu in brass sample using standard sodium thiosulfate solution (Iodometric titration)

Theory

Step 1:





Step 2:



Step 3:




Step 4:



Procedure:

0.5g of brass is dissolved in 10mL conc. HNO_3 and dilute it upto 100 mL with distilled water in standard measuring flask. Pipette 10 mL of this diluted alloy solution (test solution) into a conical flask then add a few drops of dilute NH_4OH solution until a slight permanent precipitate remains. The precipitate is redissolved by adding acetic acid drop by drop till the solution becomes clear. Then about 5 mL of a 5% KI solution is added and titrate the liberated iodine with 0.05N $\text{Na}_2\text{S}_2\text{O}_3$ solution from the burette using 1% starch solution as an indicator. The end point is marked with the disappearance of blue color. Repeat the titration 3-4 times.



OBSERVATIONS

Burette: 0.05 N $\text{Na}_2\text{S}_2\text{O}_3$ Solution

Flask: 10 ml diluted alloy solution + NH_4OH (dropwise till slight permanent Precipitation) + HAc (till the precipitation dissolves) + 5% KI




Indicator: Starch solution (2-3 drop)

End point: Disappearance of Blue color

Observation table

Sr No	Initial Burette reading (ml)	Final Burette reading(ml)	Difference (ml)	Concurrent Reading (ml)
1	0.0	9.2	9.2	8.2
2	0.0	8.2	8.2	
3	0.0	8.2	8.2	
4	0.0	8.0	8.0	

Equations:



Calculation:

$$1000 \text{ ml } 1\text{N } \text{Na}_2\text{S}_2\text{O}_3 = 63.57\text{g of Cu} = 249.57 \text{ g } \text{CuSO}_4 \cdot 5\text{H}_2\text{O}$$

$$1 \text{ ml } 1\text{N } \text{Na}_2\text{S}_2\text{O}_3 = 0.06357\text{g of Cu} = 0.24957 \text{ g } \text{CuSO}_4 \cdot 5\text{H}_2\text{O}$$

$$A \text{ ml of } 0.05 \text{ N } \text{Na}_2\text{S}_2\text{O}_3 = 8.2 \text{ ml}$$

$$\text{Quantity of Copper in the given solution} = A \times 0.05 \times 0.06357 \times 100/10$$

$$= B = 0.26 \text{ g}$$

$$\% \text{ of Copper in the alloy} = X \times 100/W = 52\%$$

Where W = weight of the alloy taken = 0.5 g

A = Constant Burette Reading and B = Quantity of Copper.

Results

1. 10 ml of diluted solution required = 8.2 ml of 0.05 N $\text{Na}_2\text{S}_2\text{O}_3$
2. Quantity of Copper in the given solution = 0.26 g
3. Percentage of Copper in the alloy = 52%.